moderna	Method Validation Report			
TITLE				
Method Validation Report: SOP-1032 Identity Confirmation of mRNA in a Lipid Nanoparticle by Sequencing Analysis				
mRNA-1273 LNP, mRNA-1273 DP				

## 1. Introduction

This report presents the method validation results of test method SOP-1032, Identity Confirmation of mRNA in a Lipid Nanoparticle by Sequencing Analysis, for mRNA-1273 LNP and DP. SOP-1032 is categorized as an identity test and results are reported as conforms or does not conform. The validation was performed at the Moderna Quality Control Laboratory following method validation protocol QC-MVP-0018 v1.0 and validation protocol addendum QC-MVP-0021 v1.0 and is based on the ICH Q2(R1) Validation of Analytical Procedures: Text and Methodology.

**SOP-1032**, Identity Confirmation of mRNA in a Lipid Nanoparticle by Sequencing Analysis, is a method used to confirm the confirmation of the open reading frame (ORF) of the mRNA LNP and DP by using Reverse transcription polymerase chain reaction (RT-PCR) to create complementary DNA template (cDNA) and analysis of cDNA by Sanger sequencing.

Method SOP-1032 was validated according to QC-MVP-0018 v1.0, using one lot of mRNA-1273 LNP and two lots of mRNA-1273 DP. The validation parameter of specify met the acceptance criteria listed in QC-MVP-0018 v1.0.

## 2. Referenced Documents

Kelefelicea Documen				
Document #	Title			
ICH Q2(R1)	International Council for Harmonization, Validation of Analytical Procedures			
FRM-0043	Read and Understand Training Form			
FRM-0742	Assay Performance Worksheet- SOP-1032 Identity Confirmation of mRNA in Drug Product by Sequencing			
SOP-1032	Identity Confirmation of mRNA in a Lipid Nanoparticle Drug Product by Sequencing Analysis			
TR-9287	- R+U - SOP-1032 and QC-MVP-0018			
TR-10486	R + U SOP-1032 QC-MVP-0021 EXS			
TR-10504	- R+U QC-MVP-0021			
QC-OTH-0199	QC-MVP-0021 Discrepancy 1			
QC-MVP-0018	Method Validation Protocol SOP-1032 ID Confirmation of mRNA in a LNP by Sequencing			
QC-MVP-0021 Method Validation Protocol SOP-1032 Identity Confirm of mRNA in LNP by Seq Analysis - Addendum				

moderna	Method Validation Report				
TITLE					
Method Validation Report: SOP-1032 Identity Confirmation of mRNA in a Lipid Nanoparticle by Sequencing Analysis					
mRNA-1273 LNP, mRNA-1273 DP					

# 3. Responsibilities

Department/ Functional Area	Responsibility
	Authors, reviews and approves validation protocols and reports.
Quality Control	Executes, reviews and approves executed data packages and data
	summaries.
	Authors validation summary reports.
O 171 A	Reviews and approves validation protocols, data summaries, and reports.
Quality Assurance	Ensures that validation documents are in alignment with Moderna
	policies and regulatory requirements.

## 4. Documentation

- 3.1. All documentation, execution, and review of the work performed for this study was conducted under current Good Manufacturing Practices (cGMP) as required by Moderna standard operating procedures.
- 3.2. Draft SOP-1032 v0.2 was followed for this testing. Assay information was documented on draft FRM-0742.
- 3.3. QC Analysts documented read and understand training on SOP-1032 v0.2 and validation protocol QC-MVP-0018 v1.0 and QC-MVP-0021 v1.0 on FRM-0043 prior to executing validation document TR-9287, TR-10486, and TR-10504 for training record. testing. Refer to
- 3.4. All relevant data collected during validation and formulae used for calculating validation characteristics will be included and discussed in this report.

# 5. Materials and Equipment

Draft SOP-1032 v0.2 and draft FRM-0742 v0.3 was followed for the validation testing.

# 5.1. Test Articles

Sample Description	Lot/Batch	RNA Concentration (mg/mL)	Summary of Analysis Document	Expected Length of RT- PCR Product (bp)
mRNA-1273 LNP <sup>1</sup>	5006820001		COA-0447	
mRNA-1273 DP <sup>2</sup>	6006820001		COA-0448	
mRNA-1273 DP <sup>2</sup>	6006920001		COA-0449	

mRNA-1273 Lipid Nanoparticle (LNP) is prepared as a liquid solution with a target mRNA concentration of

mM sodium acetate, 20 mM Tris, 87 g/L Sucrose, pH 7.5 buffer.

<sup>2</sup>mRNA-1273 Drug Product (DP) is prepared as a liquid solution with a target mRNA concentration of approximately

mg/mL in mM sodium acetate, 20 mM Tris, 87 g/L sucrose, pH 7.5 buffer.

CONFIDENTIAL - DO NOT DISTRIBUTE OUTSIDE OF MODERNATX, INC. WITHOUT WRITTEN PERMISSION.

moderna	Method Va	alidation	Report		((a)		(	(a	<b>(IIII)</b>
Method Validation Report: SOP-1032 Identity Confirmation of mRNA in a Lipid Nanoparticle by Sequencing Analysis									
mRNA-1273 LNP, mR	NA-1273 DP								

5.2. Reference Sequence of the one of the ORF of mRNA-1273

5.3. Test Primers for the of the ORF of mRNA-1273

Primer Use	Primer Name	5'-Sequence-3'
RT-PCR		
RT-PCR/ Sequencing		
Sequencing		
Sequencing		

# 5.4. Materials and Equipment

Refer to the Materials and Equipment Section of SOP-1032 v0.2.

CONFIDENTIAL - DO NOT DISTRIBUTE OUTSIDE OF MODERNATX, INC. WITHOUT WRITTEN PERMISSION.

modern	Method Validation Rep	ort	
	n Report: SOP-1032 Identity Co Sequencing Analysis	nfirmation of mRNA i	n a Lipid
mRNA-1273 LNP	mRNA-1273 DP		Shere
6.1. Validation Sum	mary Parameters and Acceptance Cri	teria	Of Valiation.
Parameter	Acceptance Criteria	Results	Pass/Fail

# 6. Validation Summary

# 6.1. Validation Parameters and Acceptance Criteria

Parameter	Acceptance Criteria	Results Kerneigh	Pass/Fail
Specificity (QC-MVP-0018)	Sequence matches 100% description of the the open reading frame with a	The sequence generated for each of three lots of mRNA-1273 matches 100% description of the of the open reading frame with	Pass
Specificity (QC-MVP-0021)	The mRNA-1273 sample sequencing matches 100% description of the the open reading frame with a the open reactions with mRNA-1273 primers does not generate sequence that 100% match the descriptions of the of the open reading frame for mRNA-1273.	The sequence generated for mRNA-1273 lot 5006820001 matches 100% description of the of the open reading frame with a  The sequence generated from non-target control with mRNA-1273 primers did not generate sequence that was a 100% match the description of the of the open reading frame for mRNA-1273.	Pass

	nderng Method Validation Report
ennice tennice money	JUETTIC
1	od Validation Report: SOP-1032 Identity Confirmation of mRNA in a Lipid particle by Sequencing Analysis
mRN	A-1273 LNP, mRNA-1273 DP
-	idation Results Specificity QC-MVP-0018 v1.0
	Specificity is the ability to establish that the test method is specific to the product-related components.
	Experimental Design:
	To assess specificity, one lot of mRNA-1273 LNP (5006820001) and two lots of mRNA-1273 DP (6006820001 and 6006920001) were prepared and analyzed per SOP-1032 v0.2 following a process.  Acceptance Criteria:
	Sequence matches 100% description of the open reading frame with a
	Results:  The sequence generated for each of three lots of mRNA-1273 matches 100% description of the
7.2.	Specificity QC-MVP-0021 v1.0
	Specificity is the ability to establish that the test method is specific to the product-related components.  Experimental Design:
	To assess specificity, one lot of mRNA-1273 LNP lot 5006820001 and one lot of non-target mRNA control lot MTDS18026 were prepared and analyzed per SOP-1032 v0.2 using primer designed for mRNA-1273 following the process.
S	Acceptance Criteria:
rent	The mRNA-1273 sample sequencing matches 100% description of the reading frame with a
	The non-target control reactions with mRNA-1273 primers does not generate sequence that 100% match the descriptions of the non-target control reactions with mRNA-1273 primers does not generate of the open reading frame for mRNA-1273.

CONFIDENTIAL - DO NOT DISTRIBUTE OUTSIDE OF MODERNATX, INC. WITHOUT WRITTEN PERMISSION.

moderna	Method Validation Report				
TITLE					
Method Validation Report: SOP-1032 Identity Confirmation of mRNA in a Lipid Nanoparticle by Sequencing Analysis					
mRNA-1273 LNP, mRNA-1273 DP					

## Results:

The sequence generated for mRNA-1273 LNP lot 5006820001 matched 100% description of the of the open reading frame with a

The sequence generated from non-target mRNA control with mRNA-1273 primers did not generate sequence that was a 100% match the description of the open reading frame for mRNA-1273.

## 8. Discrepancies

There was one discrepancy associated with the execution of this validation (QC-OTH-0199 v1.0). This discrepancy was linked to operator error and the assay run was repeated successfully. This discrepancy has no impact on this validation.

#### 9. Conclusions

Analytical test method SOP-1032 passes the acceptance criteria for all validation parameter outlined in protocols QC-MVP-0018 v1.0 and QC-MVP-0021 v1.0 for Specificity.

Analytical test method **SOP-1032** is considered validation for testing mRNA-1273 LNP and DP. A verified data summary for the validation assays is attached (**Attachment 1**).

## 10. Attachments

Attachment 1: QC-MVR-0018 Verified Data Summary

Attachment 2: ARN-20-00333-012 for QC-MVP-0018

Attachment 3: ANR-20-0401-035 for QC-MVP-0021 Invalid (

Attachment 4: ARN-20-00419-025 for QC-MVP-0021 Repeated (

# 11. Revision History

Revision #	Effective Date	Change Details	Author
1.0	Refer to Header for Effective Date	New Document	

